Laboratory Profiles

Updated 2 May, 2024

Toronto Serology Lab (TML) HbsAg, HbcAb, HCV, HIV, HTLV, CMV, Syphilis, EBV, WNV, Toxoplasmosis

Mt. Sinai Hospital 600 University Avenue, Room 1470

Toronto, ON M5G 1X5 Tel: 416 586 4800 ext. 4182

Pager: 416 380 5628 Fax: 416 619 5522

<u>Delivery Instructions:</u>

<u>Regular Hours</u> (Mon–Fri 8am to 8pm, Weekends 8am-4pm) Deliver to Microbiology Specimen Receiving Area on 14th floor (accessed by main elevators)

After Hours (Mon –Fri 8pm to 8am, Weekends 4pm-8am)
Enter using Murray Street Entrance. Sign in with security and proceed to the Murray Street elevators). Proceed to the 14th floor. Leave bags in the blue bucket in 14th floor Microbiology elevator area and ring doorbell.

Serology, NAT and West Nile testing specimen tubes required:

Donor Hospital Regions for which testing is done: Toronto, Sudbury,

• 3 clotted serum/red tubes (20 ml)

Sault Ste. Marie, Hamilton, Thunder Bay

• 3 EDTA tubes (12 ml)

TOXO Testing – done at this lab Mon-Fri 0800-1600. E-mail request to complete TOXO testing and completed requisition form to: Microbiology.SpecialQueries@sinaihealth.ca

NAT Testing **ONLY**:

• 1-2 EDTA

If NAT testing unavailable at TML, see Winnipeg lab profile below.

Out-of-Province Donor Testing:

HTLV 1/2

• 1 clotted serum/red tubes (10 ml)

West Nile Virus **ONLY**

• 1 EDTA tube (6 ml)

Notify lab/TML technologist on call one hour prior to sample arriving.

Toronto HLA Lab

HLA Laboratory - LMP University Health Network 200 Elizabeth St., <u>11E-444</u>* Toronto, ON M5G 2C4

Regular Hours (Mon-Fri 8am to 5pm)

Tel: 416 340 4995

After Hours (5pm-8am, including weekends and stat holidays)

1st call pager: 416 715 3835 2nd call pager: 416 790 4505

Fax: 416 340 3133

Delivery Instructions:

Regular Hours (Mon-Fri 8am to 5pm)

Deliver to above address

After Hours (5pm-8am, including weekends and stat holidays)

Deliver to:

HLA Laboratory – LMP °/_o Specimen Management University Health Network 200 Elizabeth St., 3E-347* Toronto, ON M5G 2C4 Transplant Programs for which testing is done: UHN, HSC, SMH

Specimen tubes required:

• Typing: 2 EDTA tubes (12 ml)

• Cross matches: 6 ACD tubes (80 ml)

Pediatric Requirements - confirm with HLA lab if needed

Call/Page HLA technologist one hour prior to sample arriving

Additional Delivery Instructions:

*Use the Elizabeth St. entrance to access the daytime and after hours drop off locations. Both drop off locations are in the Eaton Wing, accessible via the Eaton Elevators. The main door is locked between 2300-0500hrs. During this time, couriers can obtain access from security using the HELP communication box, adjacent to the wheelchair access push panel to the right of the main entrance door.

Hamilton HLA Lab

St. Joseph's Healthcare 50 Charlton Ave. East L8N 4A6

Tel: 905 522 4941 ext. 34069 Paging: 905 521 2100 ext. 76443

Delivery Instructions

Use Emergency entrance on James Street.

Deliver to Emergency Desk and have them call McMaster paging 905-521-2100 ext. 76443 - ask to have the HLA tech paged.

Transplant Program for which testing is done: St. Joe's Hamilton

Specimen tubes required: Typing

- 2 EDTA tubes (12 ml) Crossmatches
- 6 ACD tubes (80 ml)

OTDC will send blood to lab. CSC to page Tissue Typing/HLA tech on call through locating to initiate HLA

London Serology Lab (Core Lab) HbsAg, HbcAb, HCV, HIV, HTLV, Syphilis, CMV

London Health Sciences Centre, University Campus, Room A3-339 Windermere Road

Tel: 519 685 8300 ext. 35787

Fax: 519 663 3798

Core Lab to notify PRC that samples have arrived. The Core lab staff will direct the rest of the samples to the appropriate site. PRC to contact other labs as needed.

Blood from donors outside London region - sent to Core Lab University campus

London Virology Lab Stat WNV and Non-Stat EBV & Toxo

London Health Sciences Centre Victoria Hospital, Room B10-109 800 Commissioners Road East London, ON N6A 4G5

Tel: 519 685 8500 ext. 64667

Fax: 519 646 6002

1st call pager: 519 685 8500 ext. 14733

Note: EBV is tested at the Endocrinology Lab at LHSC-UH (ext.

77676)

Donor Hospital Regions for which testing is done: London, Windsor

Specimen tubes required:

- 2 SST/gold top tubes (10 ml)
- 1 EDTA tube (6 ml)
- 2 EDTA tubes (12 ml)
- 1 clotted serum/red tube (10 ml)

PRC and OTDC's – please indicate the following on specimens: After hours, couriers will leave samples with security when a security guard is on duty (20:00-06:00). Between 06:00-20:00, couriers should call the laboratory (ext. 35787) and laboratory staff will arrive immediately to pick specimens up. If at any time security is not there, please page the lab to come and

receive the sample.

Donor Hospital Regions for which testing is done: London, Windsor

Specimen tubes required:

- 1 SST/gold top tube (5 ml)
- November 1 April 30: 1 EDTA tube (6 ml)
- May 1 October 31: 2 EDTA tubes (12 ml)

Core lab to notify PRC that samples have arrived.

Nov 1 – April 30 - EBV and Toxo will be run next business day – do not page to advise samples delivered.

London HLA Lab

London Health Sciences Centre, University Campus, Room BLL 250 339 Windermere Road London, ON N6A 5A5

Tel: 519 663 3320 Fax: 519 663 3687

1st call pager: 519 685 8300 ext. 13321

Deliver samples to the SECURITY DESK at Entrance A

Donor Hospital Regions for which testing is performed: London, Windsor, Owen Sound, Sault Ste. Marie and Sudbury

Specimen tubes required: Typing

- 2 EDTA tubes (12 ml) Crossmatches
- 6 ACD tubes (80 ml)
- OTDC will notify PRC that samples have been sent to the Security Desk at UH.
- -Security staff will deliver samples to the CORE lab
- -CORE lab will call the PRC when they receive the samples

London Region Donors:

PRC will page the Primary HLA Technologist on-call toad initiate donor typing

Ottawa Serology - Children's Hospital of Eastern

Ontario (CHEO) Biochemistry Lab
HbsAg, anti-HBc total, anti-HCV, anti-HIV I/II, HIV p24,
Anti-HTVL I/II, anti-CMV, Syphilis, EBV and
Toxoplasmosis

401 Smyth Road, 3rd Floor Room 3030 Ottawa, ON K1H 8L1

Tel: 613 737 7600 ext. 2061 (alternate extension is 2069)

Fax: 613 738 4248

Donor Hospital Regions for which testing is done: Ottawa, Kingston, Belleville, Trenton

Specimen tubes required:

• 2 full clotted serum/red tubes OTDC/PRC to notify tech when specimen is on route. Send with TGLN lab requisition for stat serologies and Toxo and EBV

Drop off location - Stat Specimen Receiving Area.

Ottawa Microbiology Lab – WNV

General Campus 501 Smyth Road (3rd floor) Ottawa, ON K1H 8L6

1st Call –WNV tech on call via pager 613-402-2610 2nd call – If unable to reach tech by pager call via locating 613 737 8222

Regular hours (Mon-Fri 08:00 – 16:00)

Tel: 613 737 8322

After hours (evenings/weekends/statutory holidays)

Tel: 613 593 4135 Fax: 613 737 8324 Donor Hospital Regions for which testing is done: Ottawa, Kingston, Belleville, Trenton

Specimen tubes required:

Requirements for Hep B, Hep C, HIV Testing & West Nile

• 4 EDTA Blood Tubes

Requirements for Hep B, Hep C and HIV testing:

• 2 EDTA blood tubes

Requirements for West Nile (only completed May 1- Oct 31)

• 2 EDTA blood tubes

OTDC/PRC to page WNV tech on call for WNV testing. PRC to notify tech for non-Ottawa donors when specimen is en route

Drop off location - Stat Specimen Receiving Area The courier will have to press the buzzer to be let into the receiving area and this is where the specimens should be dropped off.

Ottawa HLA Lab Tissue Typing DNA Lab 501 Smyth Road 4th Floor, Critical Care Wing, Room 4155 Ottawa, ON K1H 8K6 Tel: 613 737 8277 Regular hours (Mon-Fri 08:00 – 16:00) Call lab directly at 613 737 8277. If no answer leave message and someone will call back. After hours (evenings/weekends/statutory holidays) Contact lab tech on-call via Ottawa Hospital locating 613 737 8899 ext. 78222	Transplant Programs for which testing is done: OHI, TOH Specimen tubes required: Typing • 2 EDTA tubes (12 ml) Crossmatches • 6 ACD tubes (48 ml) CSC to call/page Tissue Typing Tech on call if kidney accepted by TOH
Kingston HLA Lab	Transplant Program Region for which testing is done: Kingston
Kingston Immunology Lab 76 Stuart Street Located on Douglas One Monday-Friday (0800-1600): 613 549 6666 ext. 4602 After Hours (after 1600), Weekends and Holidays: 613-331-0653	Specimen tubes required: Typing • 6 ACD tubes (80 ml) Notify HLA tech to complete HLA testing and initial cross matching.
Winnipeg NAT testing Cadham Provincial Laboratory 750 William Avenue Winnipeg, MB R3E 3J7	When to send it: When the NAT testing machine at TML is out of service and NAT testing is required urgently to proceed with donation. Specimen tubes required: • 5 EDTA tubes (30 ml)
Regular hours (Mon-Fri 09:00 – 17:00 EST) Lab Tel: 204 945 0643 After hours (evenings/weekend/statutory holidays): 204 945 0643 Ask security to page the technician on-call	What to say: Please provide the rationale for NAT testing (i.e. The NAT testing machine at our Ontario lab is currently down and we have an increased risk donor that requires NAT testing results to proceed with donation)

with donation).